

Uphill Tablet Deduster C&C200G



Uphill Tablet Deduster

The machine can make production more safe and standard.

C&C200G series Uphill Tablet Deduster is new model with advanced structure. It has achieved double functions of tablet dedusting and elevating. Inlet of the machine can be connected with any model of tablet press, and outlet with metal detector, so that to realize the production line of tablet dedusting, elevating, and detecting.



Features:

- Body and cover of the machine take building-block design. Height of outlet is easy to adjust;
- 360° circle adjustment of inlet and outlet to make production more convenient and save factory room;
- Fast-joint system to ensure safe and quick disassembly and assembly;
- Outstanding efficiency can make mass production of $\Phi 3$ - $\Phi 35$ mm all-kind-shape tablet;
- Parts in contact with medicine all made of 316L tainless steel. No blind spot during whole machine cleaning. Meet cGMP demands;
- New design motor driver to make machine run more smoothly with lower noise.



C&C200G1



C&C200G1+Metal Detector

Technical Specifications	C&C200G	C&C200G1	C&C200G2	C&C200G3
Tablet Diameter (mm)	Ø3-Ø35			
Max Output Capacity (t/h)	1800000(Ø6)			
Inlet height (mm)	740-940	740-940	740-940	740-940
Outlet height (mm)	1125-1325	1525-1725	1925-2125	2325-2525
Dedusting Distance (m)	6.8	11.1	15.4	19.7
Power Supply	1PH 220V/110V 50Hz/60Hz 0.4kw			
Compression Air	0.1 m ³ /min 0.05MPa			
Vacuum	2.7 m ³ /min -0.1MPa			
Dimension (mm)	577 x 606 x (1362-1562)	577 x 606 x (1762-1962)	577 x 606 x (2162-2362)	577 x 606 x (2562-2762)
Weight (kg)	120	130	140	150

SaintyCo also offer:

Encapsulation Solutions

Tablet Compression Solutions

Coating Systems

Cartoning Systems

Blister Packing Solutions

