

GST

Gelatin Service Tank



Description

Barrels designing or manufacture is in accordance with GMP requirements of health standards, according to GB9684-88 "stainless steel food containers health standards", according to Q/WEJ01-2002 the gelatin capsule of medicine tank(container, barrel) etc. technical standards. It use three layers, which is inner, heating jacket, insulation jacket, as a whole refining by SUS304 stainless steel. Using the electric heating, intelligent temperature controller. Pressure gauge, regulator, security valve, sight glass, etc. were installed on the top of barrel. A drain valve was set in the bottom. Lid locking is the way of quick lamp. The barrel has advantages of thermal insulation, safe and convenient operation, easy cleaning and so on. This barrel has some advantages as thermal insulation, safe and convenient operation, easy cleaning and so on.

Technical Parameters	GST-100	GST-140	GST-200	GST-300
Volume(L)	100	140	200	300
Temperature	Less than 100	Less than 100	Less than 100	Less than 100
Power supply	1P 220	1P 220	1P 220	1P 220
Power(kw)	1.5	1.5	2	2.5
Pressure(Mpa)	Less than 0.05	Less than 0.05	Less than 0.05	Less than 0.05
Dimension	550*1210	600*1330	700*1330	800*1370