

Auxiliary Support Equipment

SP-200 Softgel Pulverizer



Softgel Pulverizer

Softgel Pulverizer is our patented machine. It can recycle misshaped or unqualified capsules separately by capsule shell and inside medicine solution which saves a lot cost during encapsulation process. Grinded capsule shell discharged from upper port while medicine solution comes out from lower port to avoid cross contamination. Automatic separation and discharge saves labor. It is easy clean, easy operation, easy maintenance and with simple compact structure. A desirable and efficient machine for capsule production.



SP-200

Model	SP-200
Capacity	60 kgs/hour
Grind speed	960 RPM
Net weight	80 kgs
Power supply	220V/60HZ 1.1KW
Dimension (LxWxH mm)	Φ510x1200

SaintyCo also offer:

Encapsulation Solutions

Tablet Compression Solutions

Coating Systems

Cartoning Systems

Blister Packing Solutions

